



ENGINEERING CHEMISTRY

Department of Science and Humanities

ENGINEERING CHEMISTRY

Electrochemical equilibria



Class content:

- *Glass electrode*
 - *Construction*
 - *Working*
 - *Determination of pH*

ENGINEERING CHEMISTRY

Electrochemical equilibria



Ion selective electrodes (ISE)

- Selectively respond to a **specific ion** in a mixture
- Potential developed is a function of concentration of that ion
- Have a membrane which is capable of exchanging the specific ion with solution with which it is in contact
- **Membrane electrodes**
e.g., glass electrode

ENGINEERING CHEMISTRY

Electrochemical equilibria




Glass electrode

- Ion-selective electrode
- Responds to Hydrogen ion
- pH sensitive; can be used to determine pH of a solution
- Consists of a glass membrane which is capable of exchanging H^+ ions

ENGINEERING CHEMISTRY

Electrochemical equilibria

- **Construction:**
- Glass tube , the end of which is a bulb of **very thin glass membrane**
- Glass bulb is made up of special type of glass, CORNING 015
- The glass bulb is filled with **solution of known pH** which is the reference solution
- A **silver - silver chloride electrode** is dipped inside the reference solution serves as internal reference electrode and also provides external electrical contact
- The electrode is immersed in a **solution containing H^+** which is the analyte
- **Ag/AgCl/HCl/glass**




Source: <https://glossary.periodni.com/glossary.php?en=glass+electrode>

ENGINEERING CHEMISTRY

Electrochemical equilibria

Working:

analyte solution




reference solution



Ag-AgCl electrode



- The inner and outer surfaces of the glass membrane can exchange H^+ ions with the solution they are in contact with




Source: http://www.metrohmsiam.com/teachingresearch/TRL_25/TRL25_95520780155013.pdf

- The hydrated glass membrane brings about ion exchange reaction between singly charged cations in the interstices of glass lattice and protons from the solution
- A potential is developed, which is a function of H^+ of the solution

ENGINEERING CHEMISTRY

Electrochemical equilibria

Electrode potential of a glass electrode :



boundary potential is $E_b = \frac{2.303RT}{nF} \log\left(\frac{C_1}{C_2}\right)$

since concentration of reference solution, C_2 is constant

$$E_b = L' + \frac{2.303RT}{nF} \log C_1 \quad \text{where} \quad L' = - \frac{2.303RT}{nF} \log C_2$$

At 298K, $E_b = L' + \frac{0.0591}{n} \log[H^+]$ since for H^+ , $n = 1$

$$E_b = L' - 0.0591 \text{pH}$$

ENGINEERING CHEMISTRY

Electrochemical equilibria



The glass electrode potential has 3 components

1. The boundary potential
2. The potential of internal reference electrode
3. Asymmetric potential

$$E_G = E_b + E_{ref} + E_{asymmetric}$$

- Asymmetric potential arises due to difference in responses of inner and outer surfaces of the glass bulb, due to differing conditions of stress on two glass surfaces

$$\begin{aligned} E_G &= E_b + E_{ref} + E_{asymmetric} \quad ; \quad E_b = L' + \frac{0.0591}{n} \log[H^+] \\ &= L' + \frac{0.0591}{n} \log[H^+] + E_{ref} + E_{asymmetric} \\ &= E_G^0 + 0.0591 \log[H^+] \quad \text{where } E_G^0 = L' + E_{ref} + E_{asymmetric} \\ E_G &= E_G^0 - 0.0591 \text{pH} \end{aligned}$$


ENGINEERING CHEMISTRY

Electrochemical equilibria

Determination of pH using glass electrode:

- Glass electrode is combined with an external reference electrode

$\text{Hg}/\text{Hg}_2\text{Cl}_2/\text{Cl}^-//\text{analyte solution}/\text{glass}/0.1\text{N HCl}/\text{AgCl}/\text{Ag}$



Source: <https://utkarshiniedu.wordpress.com/2016/12/22/lecture-1-108-ion-selective-electrodes/>

ENGINEERING CHEMISTRY

Electrochemical equilibria

Determination of pH using glass electrode:

- The emf of the cell is determined potentiometrically

$$E_{\text{cell}} = E_{\text{G}} - E_{\text{calomel}} ; \quad E_{\text{G}} = E^0_{\text{G}} - 0.0591\text{pH}$$

$$= E^0_{\text{G}} - 0.0591\text{pH} - E_{\text{calomel}}$$

$$\text{pH} = \frac{E^0_{\text{G}} - E_{\text{calomel}} - E_{\text{cell}}}{0.0591}$$

- To evaluate E^0_{G} the glass electrode is dipped in a **solution of known pH(buffer solution)** and combined with calomel electrode, the emf of the cell is measured from which **E^0_{G} can be evaluated**

ENGINEERING CHEMISTRY

Electrochemical equilibria

Applications of glass electrode:

Used extensively in chemical, industrial, agricultural and biological labs

Advantages of glass electrode :


- Can be used in oxidizing and reducing environments and metal ions
- Does not get poisoned
- Can be used for very small volumes
- Accurate results can be obtained between pH 1 to 9 by ordinary electrodes. However by using special glass electrodes pH 1 to 14 can be measured with accuracy
- Simple to operate and can be used with portable instruments

ENGINEERING CHEMISTRY

Electrochemical equilibria


Disadvantages of glass electrode:

- Because of high resistance of glass, a simple potentiometer cannot be used. It requires sensitive potentiometer for emf measurements
- Glass membrane is very delicate, hence has to be handled carefully
- At very high pH levels usually over a pH of 9, **Alkaline error** is observed



When the Sodium ion level is relatively high, some of the H^+ ions in the gel layer around the sensitive electrode membrane are replaced by Na^+ ions

The electrode may eventually respond to Na^+ instead of H^+ ions, giving a false lower pH value than the actual value





THANK YOU

Department of Science and Humanities